

RE-V/AP/S WINE CLARIFIERS

REDA
SEPARATION

Model	Grape juice	Wine at end fermentation	Pre-decanted wine	Sparkling wine	
				< 2.5 bar	> 2.5 bar
RE30V	2,000 - 3,000	4,000 - 6,000	5,000 - 7,000	3,000 - 4,000	----
RE30S	2,000 - 3,000	4,000 - 6,000	5,000 - 7,000	3,000 - 4,000	3,000 - 4,000
RE50V	3,000 - 5,000	5,000 - 8,000	6,000 - 9,000	5,000 - 6,000	----
RE50S	3,000 - 5,000	5,000 - 8,000	6,000 - 9,000	5,000 - 6,000	5,000 - 6,000
RE85V/AP	6,000 - 10,000	12,000 - 15,000	15,000 - 20,000	10,000 - 12,000	----
RE85S	6,000 - 10,000	12,000 - 15,000	15,000 - 20,000	10,000 - 12,000	10,000 - 12,000
RE130V/AP	10,000 - 14,000	18,000 - 25,000	20,000 - 26,000	15,000 - 18,000	----
RE130S	10,000 - 14,000	18,000 - 25,000	20,000 - 26,000	15,000 - 18,000	15,000 - 18,000

- Flow rates are expressed in liters/hour
- AP version are supplied with "self-thinking" system

